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Digital Camera Image Display／
Measurement Program Installation／
Operations Manual
The 6th Edition

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・ Installation of DirectX9

- (1) "dxsetup.exe" in the DirectX9 folder of CD-ROM for installation is performed.
- (2) It chooses "to agree" as the consent contract of a screen, and the "Next" is chosen.
- (3) The "Next" is chosen and installation is started.
- (4) If completion of installation is displayed, a screen will be closed and a personal computer will be rebooted.

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・ Installation of a camera driver

- (1) "Setup.exe" in the Driver folder of CD-ROM for installation is performed.
- (2) It chooses from "Next" in order of "install", and a driver is installed.
- (3) A USB camera will be connected if installation is completed.

Connection of a camera will start installation of a driver automatically. I carry out

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・ Installation of a USB key driver

- (1) "InstDrv.exe" in the USBKEY folder in CD-ROM for installation is performed.
- (2) It chooses from "Next" in order of "Next", and a driver is installed.
- (3) "Complete" will be chosen if the completion screen of installation is displayed.
- (4) A personal computer is rebooted.

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・ Installation of a digital camera image display / measurement program

- (1) "Install.exe" in CD-ROM for installation is performed.
- (2) Since an installation screen is displayed, the "Next" is chosen.
- (3) Since the screen which chooses the folder of an installation place is displayed, the "Next" is chosen.
- (4) If installation is completed, since an end screen will be displayed, "Completion" is chosen.
- (5) Installation is completion.

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
button will end a program. In this case, keep in mind that no measurement result, image, etc. are saved.

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Keep in mind that it will end while a result has not been saved if "yes" is chosen.

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● Explanation of a function of the screen upper part

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It is movie mode at the time of starting.

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When measuring, it is made still mode.

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When changing the mode from a still picture to a movie and there is a non-saved measurement result, a movie mode change check is displayed. Keep in mind that a result is not saved if "yes" is chosen.

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
The still picture file (BMP, JPG form) saved is opened.

If a file is opened, it will change to still picture mode automatically.

The image size of a still picture file is too large to display. A check message is displayed when it cannot display.

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A still picture (BMP, JPG form) and a movie (AVI form) are saved.

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○ The method of save of a movie

"Save" is chosen and a save place is specified.

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A file name is inputted and the kind of file is made into "a movie (*.avi)."

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
If "save" is chosen, the icon of "movie" will change to "start".

"start" is pushed and save of a movie is started.

If the remaining time which chooses a "stop" for stopping or is displayed on the comment field is set to 0, it will stop automatically. (The capacity of about 2GB is consumed in 1 minute as a standard)

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
It is an option set up at the time of save of a still picture.


"JPEG quality" sets up the quality of image at the time of save.

If a number is clicked, the screen of quality-of-image adjustment will come out.

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"Result save" puts in a check to save still pictures also including the measured result.

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It is used to indicate the image currently displayed by expansion.

The range of the magnification which can be specified is from x0.5 to x2.0.

When the image is too large, it may not indicate by expansion. Image size is set to x1.0 after a confirmation screen is displayed, when not indicating by expansion.

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A setup of "brightness", "contrast", etc., etc. can be performed at the time of movie mode.

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
A comment can be freely inputted to two lines. The inputted comment is saved as a text file of a same name at the time of still picture save. Also when image is opened, it is read simultaneously and displayed. Moreover, the operation method of the measurement specified during measurement is displayed. The operation method and the remaining time are displayed at the time of save of a movie.

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It can be used by still picture mode. The range and area are measurable.

The icon of an item to measure is clicked on it and used.

The icon under operation is displayed in yellow.

A measurement result is displayed on the measurement value at the lower left of a screen.

During specification of a point or the range, if it right-clicks on a still picture, it will return at the operation before one.

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- ・ Distance(two-point directions)/ (three-point directions)

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The specified distance of two points is measurable.

- ・ Two-point directions

The starting point is specified by the 1st click and a terminal point is specified by the 2nd click.

- Three-point directions

Two places, the starting point and a terminal point, are clicked first, and a base line is drawn.

In accordance with a motion of a mouse, a base line, a parallel line, and the line that connects two lines are displayed.

A click of the 3rd point measures the distance between a base line and parallel lines. When the outside of a base line is clicked, an error message is displayed, and it returns to specification of the 3rd point.

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The diameter of the circle which passes along three specified points is measured.

Since a circle will be displayed according to cursor if the two circumferences of a circle to measure are clicked, the 3rd point in all is specified as a circle to measure.

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The distance from each diameter of two circles and the center of two circles which were specified is measured.

The 3rd point in all is specified as a circle to measure if the two circumferences of a circle to measure are clicked, since a circle will be displayed according to cursor. Specification of two circles will measure the distance of two circles automatically.

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
It will be underlined with the altitude to a base line, if two places, the starting point and a terminal point, are clicked first, a base line is drawn and a corner point is specified.

A corner point can be specified continuously.

A right-click will end measurement of Per line.

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
The area of the specified range is measured.

Since a confirmation screen will be displayed if the range to measure is specified and the same position as the starting point is specified, "yes" is pushed and it measures.

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・ Angle (three-point directions)/Angle 2 (four-point directions)

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The angle of the specified part is measured.

The measurement method specifies either of two kinds, three-point directions and four-point directions, on a setting screen.

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A line will be displayed if the 2nd point is clicked with the 1st point.

Next, specification of the 3rd point will calculate the inside.

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The first line will be displayed if the 2nd point is clicked with the 1st point.

The second line will be displayed if the 4th point is clicked with the 3rd point.

Specification of the side which measures an angle will calculate the inside or the outside corner.

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It is counted how many the point is specified and it specified.

The number counted when the count part is clicked and the right-click was finally pushed is displayed.

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
It clicks and fixes in a part to specify since cross line are displayed on a screen.

A button will be pushed in during the display of cross line.

If it clicks on an icon while on display, the display of cross line will disappear.

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
Two cross line can be specified by setup.

It clicks and fixes in a part to specify since the 1st cross line are displayed on a screen. Then, since the 2nd cross line are displayed, it clicks and fixes in a part to specify similarly.

A button will be pushed in during the display of cross line.

If it clicks on an icon while on display, the display of cross line will disappear.

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A scale is displayed on a screen.

A setting screen is displayed first. If the length of X and Y, the display style of a scale, and a display position are specified and the O.K. button is clicked, a scale will be displayed on a screen.

If it clicks on an icon while on display, the display of a scale will disappear.

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
A click of a "setting" button will display a measurement function-settings screen.

Change of a color can be performed if the color of the item which a measurement line and a measurement value want to change is clicked.

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Cross line can perform change of a color as well as selection of a number, and other lines.

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The scale can perform change of a color.

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The display below the decimal point of a measurement value can be set up from 0 to 8 figures.

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The character size of a measurement value and a scale can be set up.

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
The method (three-point directions or four-point directions) of angle measurement can be set

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The method (two-point directions or three-point directions) of measurement can be set up for two points.

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
“OK” is chosen and the contents of change are made to reflect.

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If a proofreading value is set up, measurement united with the absolute size can be performed.

A proofreading value can be set up to ten pieces.


When changing a setup of a proofreading value, the setting button beside a proofreading value is clicked. Since it is cleared when all measuring at this time, the warning screen is displayed.

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Although the absolute size is known, two places, the starting point and a terminal point, are specified.

Since the following screen will be displayed if a terminal point is specified, the length (standard length) of a name and the specified line and a unit are inputted.


A dialog box titled "Setting" with a blue header. It contains three input fields: "Name" with the text "Proof 1", "STD length" with the value "1", and "Unit" with a dropdown menu showing "MIC". At the bottom are two buttons: "OK" and "Cancel".

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If "OK" is chosen and set up, it will be measured according to the fiducial point.

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- Measurement value

A dialog box titled "Measurement value" with a yellow header. It contains three icons: a floppy disk for "Save", a red X for "Delete", and a trash can for "Clear". Below these is a table with three columns: "Type", "Value", and a unit column. The table has several rows of data, with the first row highlighted in blue.

Type	Value	
Distance	193.65	MIC
Diameter	134.14	MIC
Count	4	
Angle	53.46	DEG
Area	2751.50	MIC2
Per line	79.85	MIC

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The result calculated by "measurement" is displayed. Selection of the item of a result will change the color of the corresponding section on a still picture side.

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A measurement result can be outputted as a CSV file.

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The unnecessary contents of measurement can be deleted.

Contents to delete are chosen and a deletion icon is selected.

Multiple selection of the continuous contents of measurement can be made, and they can also be deleted.

Since a confirmation screen comes out, "yes" is chosen and it performs.

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All the measured contents are deleted.

Since a confirmation screen comes out, "yes" is chosen and it performs.